

HELINI

Y microdeletion

Real-time PCR Kit

AZF-A	AZF-B	AZF-C
sY746/sY81/ sY84/sY86	sY118/sY113/ sY127/sY134/sY142	sY153/sY157/sY158 /sY254 /sY255

Instructions for use

For use with: Agilent, Bio-Rad, Roche Lightcycler-96, Roche-Z480/Cobas-480, Applied Bio systems [ABI], Thermo-Piko-Real, Rotor gene 5/6plex, Alta-96, Cepheid Real time PCR machines.



8604



25



HELINI Biomolecules, Chennai, INDIA

Intended Use

The HELINI Y Chromosome Microdeletion Real-time PCR Kit is an in vitro nucleic acid amplification kit for the detection of 14 micro deletions in human Y chromosome [AFZ – A/B/C region] associated with male infertility.

Kit components

Components	Volume Per reaction	Volume
Probe PCR Master Mix	10µl	250µl x 15
sY81 PP Mix	5µl	125µl
sY84 PP Mix	5µl	125µl
sY86 PP Mix	5µl	125µl
sY746 PP Mix	5µl	125µl
sY118 PP Mix	5µl	125µl
sY113 PP Mix	5µl	125µl
sY127 PP Mix	5µl	125µl
sY134 PP Mix	5µl	125µl
sY142 PP Mix	5µl	125µl
sY153 PP Mix	5µl	125µl

Components	Volume Per reaction	Volume
sY157 PP Mix	5µl	125µl
sY158 PP Mix	5µl	125µl
sY254 PP Mix	5µl	125µl
sY255 PP Mix	5µl	125µl
SRY14 PP Mix	5µl	125µl
Positive control Mix	10µl	150µl
PCR grade water	1	4ml

Storage

- The kit is shipped on gel ice. Upon arrival, all components should be stored in -20°C. They are stable until the expiration date stated on the label.
- Repeated thawing and freezing should be avoided, as this might affect the performance of the assay.
- If the reagents are to be used only intermittently, they should be frozen in aliquots. Storage at 2 to 8°C should not exceed a period of 5 hours.

Material and instruments required

- Real-time PCR instrument having FAM channel
- Automatic Nucleic acid extraction system or spin column based purification kit for the purification of nucleic acids
- Desktop centrifuge having 13000rpm or above with a rotor for 1.5/2 ml reaction tubes
- Centrifuge with a rotor for PCR strips/tubes and 96 well plates
- Optical cap qPCR tubes or strips or 96 wells
- Micro Pipettes (variables)
- Micro Pipette tips with filters (disposable)
- Powder-free gloves (disposable)

[Please ensure that all instruments used have been installed, calibrated, checked and maintained according to the manufacturer's instructions and recommendations.]

Product Use Limitations

- All reagents may exclusively be used in molecular diagnosis.
- The product is to be used by personnel specially instructed and trained in Molecular diagnosis.
- Strict compliance with the user manual is required for optimal PCR results.
- Attention should be paid to expiration dates printed on the box and labels of all components. Do not use expired components.
- Wear protective disposable powder-free gloves, a laboratory coat and eye protection when handling specimens and kit components.
- Avoid microbial and nuclease (RNase/RNase) contamination of the specimens and the components of the kit.
- Always use RNase/RNase-free disposable pipette tips with aerosol barriers.
- Use separated and segregated working areas for sample preparation, reaction setup and amplification/detection activities.
- The workflow in the laboratory should proceed in unidirectional manner. Always wear disposable gloves in each area and change them before entering a different area.
- Store positive and/or potentially positive material separated from all other components of the kit.
- Do not open the reaction tubes/plates post amplification, to avoid contamination with amplicons.

- Additional controls may be tested according to guidelines or requirements of local, state and/or federal regulations or accrediting organizations.
- Do not autoclave reaction tubes after the PCR, since this will not degrade the amplified nucleic acid and will bear the risk to contaminate the laboratory area.
- Discard sample and assay waste according to your local safety regulations.

Technical Assistance

For technical assistance and more information, please contact;

0091-44-244490433

helinibiomolecules@gmail.com

Product description

HELINI Y Chromosome Microdeletion Real-time PCR Kit constitutes a ready-to-use system for the detection of microdeletion in human Y chromosome using polymerase chain reaction (PCR). It contains reagents and enzymes for the specific amplicon in fluorescence channels FAM. Human endogenous control [sY14] is provided to justify the DNA purification and PCR inhibition. External mutation positive control is supplied to assist the run.

Specificity

Y chromosome microdeletion primer and probe have been designed for the specific and exclusive *in vitro* detection of following microdeletion.

AZF-A	AZF-B	AZF-C
sY746/sY81/ sY84/sY86	sY118/sY113/ sY127/ sY134/sY142	sY153/sY157/sY158/ sY254 /sY255

The primers and probe sequences in this kit have 100% homology with clinically relevant reference sequences based on a comprehensive bioinformatics analysis.

Dynamic linear range

The linear range was evaluated by analyzing a logarithmic dilution series of DNA concentrations ranging from 100ng/μl to 10ng/μl. At least six replicates per dilution were analyzed. The slopes are in expected limit in the recommended DNA concentration of 10ng/μl.

Analytical Sensitivity

The analytical sensitivity is defined as the concentration of DNA molecules (ng/μl) that can be detected with a positivity rate of 95%. The analytical sensitivity was determined by analysis of dilution series of quantified MPL mutations specific DNA from 1ng/μl to 100ng/μl in triplicates. Under optimal PCR conditions, the analytical sensitivity is 2ng/μl

Note:

DNA Purification

Purified DNA is the starting material for the Real-time PCR assay. The quality of the purified DNA has a profound impact on the performance of the entire test system. It has to be ensured that the purification system used for DNA purification is compatible with real-time PCR technology.

If you are using a spin column based sample preparation procedure having washing buffers containing ethanol, it is highly recommended to perform an additional centrifugation step for 3min at approximately 17000 x g (~ 13000 rpm), using a new collection tube, prior to the elution of the DNA.

Detection Protocol

Things to do before starting

- Before use, all kit components need to be thawed completely, mixed by gently inverting and centrifuged briefly.
- Make sure that Positive and Negative control is included in every run.

Components	Volume per reaction
Probe PCR Master Mix	10μl
PP Mix	5μl
Purified DNA	10μl
Final reaction volume	25μl

Negative Control setup [NTC]

Add 10μl of PCR grade water.

Positive Control setup

Add 10μl of the Positive control

Centrifuge PCR vials briefly before placing into thermal cycler.
[Note: There should not be any bubbles in the reaction mix. Bubbles interfere with fluorescence detection.]

Programming Thermal cycler

Sample volume	25µl
Fluorescence Dyes	FAM
Passive reference	None
Ramping rate	Default

Thermal Profile

	Step	Time	Temp
	Taq enzyme activation / Hold	15min	95°C
35 cycles	Denaturation	20sec	95°C
	Annealing/Data collection*	20sec**	60°C
	Extension	20sec	72°C

** Some qPCR machines may require minimum 30sec for data collection; in that case, set to 30sec, this will not affect the performance.

Data collection/Acquisition	Targets
FAM	All targets & endogenous control-SRY14

Reading the graph:**Endogenous control**

Human Y specific gene [sY14] is given as endogenous control. It amplifies a single copy human gene from the test samples. A successful amplification indicates that test sample is properly collected and has its biological property with required number of human cells for PCR. The endogenous control CT value must be in range of 20 +/- 9.

Step-2 –Negative and Positive control validation

Select the NTC and Positive control and view the graph of amplification.

The NTC must be flat with no Ct value. If required adjust the threshold value just above the NTC. The Positive control must be amplified.

NTC justifies NO contamination in the reagent as well as fine pipetting and its environment. PC justifies the reagents storage conditions and reaction parameters are as prescribed.

Step-3 Test Sample status

In FAM channel, select test sample well one by one, analyze the graph/amplification.

Qualitative interpretation of results:

FAM-Amplified					
Targets	NTC	PC	Endo/ SRY14	Sample	Interpretation
AZF-A					
sY746	No	Yes	Yes	Yes	No deletion
sY81	No	Yes	Yes	Yes	No deletion
sY84	No	Yes	Yes	Yes	No deletion
sY86	No	Yes	Yes	Yes	No deletion
AZF-B					
sY118	No	Yes	Yes	Yes	No deletion
sY113	No	Yes	Yes	Yes	No deletion
sY127	No	Yes	Yes	Yes	No deletion
sY134	No	Yes	Yes	Yes	No deletion
sY142	No	Yes	Yes	Yes	No deletion
AZF-C					
sY153	No	Yes	Yes	Yes	No deletion
sY157	No	Yes	Yes	Yes	No deletion
sY158	No	Yes	Yes	Yes	No deletion
sY254	No	Yes	Yes	Yes	No deletion
sY255	No	Yes	Yes	Yes	No deletion

Note:

Absence of amplification in target indicates a micro deletion.

No: Not amplified

Yes: Amplified

Limitations

Good laboratory practice is essential for proper performance of this assay. Strict compliance with the instructions for use is required for optimal results.

Analysts should be trained and familiar with testing procedures and interpretation of results prior to performing the assay.

A false negative result may occur due to improper collection, transport or handling. Appropriate specimen collection, transport, storage and processing procedures are required for the optimal performance of this test.

Extreme care should be taken to preserve the purity of the components of the kit and reaction setups. All reagents should be closely monitored for impurity and contamination. Any suspicious reagents should be discarded.

The presence of PCR inhibitors may cause under quantification, false negative or invalid results.

Potential mutations within the target regions of the genome covered by the primers and/or probes used in the kit may result in under quantification and/or failure to detect.

As with any diagnostic test, the HELINI Y Microdeletion Real-time PCR results need to be interpreted in consideration of all clinical and laboratory findings.

Quality Control

In accordance with the HELINI Biomolecules in house Quality Management System, each lot of HELINI Y Microdeletion Real-time PCR kit is tested against predetermined specifications to ensure consistent product quality.

Explanations of symbols



In vitro diagnostic medical device



Catalogue number



Pack size – number of tests



Manufacturer

Manufactured by

HELINI Biomolecules,

Ohmlina, 26, 2nd Avenue,

Khuthubi Complex, Vettuvankeni,

Chennai 600115, Tamilnadu, INDIA

www.helini.in

info@helini.in

helinibiomolecules@gmail.com

+91-44-24490433

