

HELINI

BRAC1 & BRAC2

Real-time PCR Kit

Gene	Mutations
BRCA1	185delAG / 300T>G [C61G] / 2080delA 4153delA / 5382insC
BRCA2	6174delT

Instructions for use

For use with: Agilent, Bio-Rad, Roche Lightcycler-96, Roche-Z480/Cobas-480, Applied Bio systems [ABI], Thermo-Piko-Real, Rotor gene 5/6plex, Alta-96, Cepheid Real time PCR machines.



8510



25



HELINI Biomolecules, Chennai, INDIA

Contents

Intended use	5
Kit components	5
Storage	6
Material and Instruments required	6
Product use limitations	7
Product description	10
Specificity & Linear range	10
Analytical sensitivity	11
Note – DNA purification	11
Detection protocol – reaction setup	12
Programming – Thermal profile & Dye selection	13
Reading the graph	15
Qualitative interpretation	17
Limitations	18
Quality control	19

Intended Use

The HELINI BRCA1 & BRCA2 Mutations Real-time PCR Kit is an in vitro nucleic acid amplification kit for the detection of 6 somatic mutations in the BRCA1 & BRCA2 genes associated with high risk of breast and/or ovarian cancer.

Kit components

Components	Volume Per reaction	No of vials	Reaction per vials
185delAG Master Mix	15µl	2	17
300T Master Mix	15µl	2	17
300G Master Mix	15µl	2	17
2080delA Master Mix	15µl	2	17
4153delA Master Mix	15µl	2	17
5382insC Master Mix	15µl	2	17
6174delT Master Mix	15µl	2	17
Positive control Mix	10µl	1	250µl
Water, PCR grade		1	4ml

Storage

- The kit is shipped on gel ice. Upon arrival, all components should be stored in -20°C. They are stable until the expiration date stated on the label.
- Repeated thawing and freezing should be avoided, as this might affect the performance of the assay.
- If the reagents are to be used only intermittently, they should be frozen in aliquots. Storage at 2 to 8°C should not exceed a period of 5 hours.

Material and instruments required

- Real-time PCR instrument having FAM channel
- Automatic Nucleic acid extraction system or spin column based purification kit for the purification of nucleic acids
- Desktop centrifuge having 13000rpm or above with a rotor for 1.5/2 ml reaction tubes
- Centrifuge with a rotor for PCR strips/tubes and 96 well plates
- Optical cap qPCR tubes or strips or 96 wells
- Micro Pipettes (variables)
- Micro Pipette tips with filters (disposable)
- Powder-free gloves (disposable)

[Please ensure that all instruments used have been installed, calibrated, checked and maintained according to the manufacturer's instructions and recommendations.]

Product Use Limitations

- All reagents may exclusively be used in molecular diagnosis.
- The product is to be used by personnel specially instructed and trained in Molecular diagnosis.
- Strict compliance with the user manual is required for optimal PCR results.
- Attention should be paid to expiration dates printed on the box and labels of all components. Do not use expired components.
- Wear protective disposable powder-free gloves, a laboratory coat and eye protection when handling specimens and kit components.
- Avoid microbial and nuclease (RNAse/RNase) contamination of the specimens and the components of the kit.
- Always use RNAse/RNase-free disposable pipette tips with aerosol barriers.
- Use separated and segregated working areas for sample preparation, reaction setup and amplification/detection activities.
- The workflow in the laboratory should proceed in unidirectional manner. Always wear disposable gloves in each area and change them before entering a different area.
- Store positive and/or potentially positive material separated from all other components of the kit.
- Do not open the reaction tubes/plates post amplification, to avoid contamination with amplicons.

- Additional controls may be tested according to guidelines or requirements of local, state and/or federal regulations or accrediting organizations.
- Do not autoclave reaction tubes after the PCR, since this will not degrade the amplified nucleic acid and will bear the risk to contaminate the laboratory area.
- Discard sample and assay waste according to your local safety regulations.

Technical Assistance

For technical assistance and more information, please contact;

0091-9382810333

0091-44-244490433

helinibiomolecules@gmail.com

Product description

HELINI BRCA1 & BRCA2 Real-time PCR Kit constitutes a ready-to-use system for the detection of BRCA1 & BRCA2 gene mutations using polymerase chain reaction (PCR). It contains reagents and enzymes for the specific amplification for the direct detection of the specific amplicon in fluorescence channels FAM. External positive control is supplied to assist the run.

Specificity

BRCA1 & BRCA2 mutation primer and probe have been designed for the specific and exclusive *in vitro* detection of BRCA1 & BRCA2 mutations.

Gene	Mutations
BRCA1	185delAG / 300T>G [C61G] / 2080delA 4153delA / 5382insC
BRCA2	6174delT

The primers and probe sequences in this kit have 100% homology with clinically relevant reference sequences based on a comprehensive bioinformatics analysis.

Dynamic linear range

The linear range was evaluated by analyzing a logarithmic dilution series of mutation specific pDNA concentrations ranging from 1.00E+09 to 1.00E+00copies/μl. At least six replicates per dilution were analyzed. The linear range is 1.00E+06 to 1.00E+00copies/μl.

Analytical Sensitivity

The analytical sensitivity is defined as the concentration of plasmid DNA molecules (copies/μl) that can be detected with a positivity rate of 95%. The analytical sensitivity was determined by analysis of dilution series of quantified mutation specific pDNA from 0.001copies to 10copies/μl in triplicates. Under optimal PCR conditions, the analytical sensitivity is 1 copy per micro liter.

Note:**DNA Purification**

Purified DNA is the starting material for the Real-time PCR assay. The quality of the purified DNA has a profound impact on the performance of the entire test system. It has to be ensured that the purification system used for DNA purification is compatible with real-time PCR technology.

If you are using a spin column based sample preparation procedure having washing buffers containing ethanol, it is highly recommended to perform an additional centrifugation step for 5 min at approximately 17000 x g (~ 13000 rpm), using a new collection tube, prior to the elution of the DNA.

Detection Protocol

Things to do before starting

- Before use, all kit components need to be thawed completely, mixed by gently inverting and centrifuged briefly.
- Make sure that Positive and Negative control is included in every run.
- Include 0.5 reaction volume for pipetting error while calculating the volume for total number of reactions.

In a sterile 0.1 or 0.2ml qPCR tube, add								
Components	185	300T	300G	2080	4153	5382	6174	Endo
Master Mix	15µl	15µl	15µl	15µl	15µl	15µl	15µl	15µl
Mix DNA with master mix by pipetting up and down								
DNA [>5ng/µl]	10µl	10µl	10µl	10µl	10µl	10µl	10µl	10µl
Final rxns volume	25µl	25µl	25µl	25µl	25µl	25µl	25µl	25µl

Negative Control setup [NTC]

Add 10µl of PCR grade water.

Positive Control setup

Add 10µl of the Positive control

Centrifuge PCR vials briefly before placing into thermal cycler.

[Note: There should not be any bubbles in the reaction mix. Bubbles interfere with fluorescence detection.]

Programming Thermal cycler

Sample volume	25µl
Fluorescence Dyes	FAM
Passive reference	None
Ramping rate	Default

Thermal Profile

	Step	Time	Temp
	Taq enzyme activation / Hold	15min	95°C
35 cycles	Denaturation	20sec	95°C
	Annealing/Data collection*	20sec	60°C
	Extension	20sec	72°C

Data collection/Acquisition	Targets
FAM	185delAG / 300T>G [C61G] / 2080delA 4153delA / 5382insC
	BRCA2 - 6174

Reading the graph:**Step-1 –Negative and Positive control validation**

Select the NTC and Positive control and view the graph of amplification.

The NTC must be flat with no Ct value. If required adjust the threshold value just above the NTC. The Positive control must be amplified.

NTC justifies NO contamination in the reagent as well as fine pipetting and its environment. PC justifies the reagents storage conditions and reaction parameters are as prescribed.

Step-2 Test Sample status

In FAM channel, select test sample well one by one, analyze the graph/amplification.

Qualitative interpretation of results:

Gene	Mutations	NTC	Positive control	Sample amplified	Sample Not Amplified
BRCA1	185delAG	Negative	Positive	No mutation Detected	Mutation Detected
BRCA1	300T	Negative	Positive	Wild	Mutation
BRCA1	300G	Negative	Positive	Mutation	Wild
BRCA1	2080delA	Negative	Positive	No Mutation Detected	Mutation
BRCA1	4153delA	Negative	Positive	No Mutation Detected	Mutation
BRCA1	5382insC	Negative	Positive	No Mutation Detected	Mutation
BRCA2	6174delT	Negative	Positive	No Mutation Detected	Mutation

Note: Interpretation for BRCA1-300T>G

300T	300M	Sample Status
<i>Positive</i>	Negative	Homozygous T/T [Wild]
Negative	<i>Positive</i>	Homozygous G/G [Mutant]
<i>Positive</i>	<i>Positive</i>	Heterozygous T/G

Limitations

Good laboratory practice is essential for proper performance of this assay. Strict compliance with the instructions for use is required for optimal results.

Analysts should be trained and familiar with testing procedures and interpretation of results prior to performing the assay.

A false negative result may occur due to improper collection, transport or handling. Appropriate specimen collection, transport, storage and processing procedures are required for the optimal performance of this test.

Extreme care should be taken to preserve the purity of the components of the kit and reaction setups. All reagents should be closely monitored for impurity and contamination. Any suspicious reagents should be discarded.

The presence of PCR inhibitors may cause under quantification, false negative or invalid results.

Potential mutations within the target regions of the genome covered by the primers and/or probes used in the kit may result in under quantification and/or failure to detect.

As with any diagnostic test, the HELINI BRCA1 & BRCA2 Real-time PCR results need to be interpreted in consideration of all clinical and laboratory findings.

Quality Control

In accordance with the HELINI Biomolecules in house Quality Management System, each lot of HELINI BRCA1 & BRCA2 Real-time PCR kit is tested against predetermined specifications to ensure consistent product quality.

Explanations of symbols



In vitro diagnostic medical device



Catalogue number



Pack size – number of tests



Manufacturer

Manufactured by

HELINI Biomolecules,
Ohmlina, 26, 2nd Avenue,
Khuthubi Complex, Vettuvankeni,
Chennai 600115, Tamilnadu, INDIA
www.helini.in
info@helini.in
helinibiomolecules@gmail.com
+91-44-24490433
+91-9382810333

