

# HELINI

## Purefast

### Stool processing buffer

Instructions for use

**For use with:** Human stool samples

CE

IVD

REF

2009



100ml



HELINI Biomolecules, Chennai, INDIA

## Contents

Intended use	5
Kit components	5
Storage	6
Material and Instruments required	6
Purification protocol	9
Quality control	10

**Intended Use**

HELINI Stool processing buffer is designed for rapid and cost-effective small-scale preparation of high-quality DNA/RNA from Stool samples. It preserves nucleic acids as well as removes inhibitor from stool samples. Purified DNA/RNA can be used directly in PCR/qPCR and other molecular biology enzymatic reactions.

**Kit components**

Components	Volume Per reaction	Volume
Stool processing buffer	2ml	100ml

**Storage**

- The kit is shipped in room temperature.
- Store at room temperature.
- They are stable until the expiration date stated on the label.
- Repeated thawing and freezing should be avoided, as this might affect the performance of the assay.

**Material and instruments required**

- 5ml micro centrifuge tubes
- Desktop centrifuge having 13000rpm or above with a rotor for 1.5/2/5ml reaction tubes
- Micro Pipettes (variables)
- Micro Pipette tips with filters (disposable)
- Powder-free gloves (disposable)

*[Please ensure that all instruments used have been installed, calibrated, checked and maintained according to the manufacturer's instructions and recommendations.]*

**Product Use Limitations**

- All reagents may exclusively be used in molecular biology DNA/RNA applications.
- The product is to be used by personnel specially instructed and trained in Molecular biology experiments.
- Strict compliance with the user manual is required for optimal PCR results.
- Attention should be paid to expiration dates printed on the box and labels of all components. Do not use expired components.
- Wear protective disposable powder-free gloves, a laboratory coat and eye protection when handling specimens and kit components.
- Avoid microbial and nuclease (DNase/RNase) contamination of the specimens and the components of the kit.
- Always use DNase/RNase-free disposable pipette tips with aerosol barriers.
- Use separated and segregated working areas for sample preparation, reaction setup and amplification/detection activities.
- The workflow in the laboratory should proceed in unidirectional manner. Always wear disposable gloves in each area and change them before entering a different area.
- Store positive and/or potentially positive material separated from all other components of the kit.
- Discard sample and assay waste according to your local safety regulations.

**Technical Assistance**

For technical assistance and more information, please contact;  
0091-9382810333  
0091-44-244490433  
helinibiomolecules@gmail.com

**Procedure:****For storage:**

1. Transfer 2ml of Stool Processing buffer into 5ml centrifuge tube.
2. Transfer 50 to 150mg of stool samples and Vortex thoroughly for 2 to 5min. [Make sure that stool sample mixed thoroughly in Stool processing buffer]
3. Incubate at room temperature for 5min and Store at -20C for future use. For immediate use, please continue.

**For purification**

4. Mix well by inverting several times and Centrifuge at 13000rpm for 5min
5. Transfer 350µl supernatant into fresh 2ml centrifuge tube. [Note: There may be floating particles in the top of the supernatant. Do not collect floating particles. Insert micro tip under floating layer and collect the clear supernatant]
6. Use this 350µl of supernatant for the DNA/RNA purification. Complete the purification process as per the kit manufacture instructions.

Note: Frozen sample has to be thawed to room temperature and vortexed thoroughly for 5 min and follow steps from 3.

**Quality Control**

In accordance with the HELINI Biomolecules in house Quality Management System, each lot of HELINI Purefast Stool processing buffer is tested against predetermined specifications to ensure consistent product quality.

**Explanations of symbols**

In vitro diagnostic medical device



Catalogue number



Pack size – number of tests



Manufacturer

Manufactured by

**HELINI *Biomolecules***,

Ohmlina, 26, 2<sup>nd</sup> Avenue,

Khuthubi Complex, Vettuvankeni,

Chennai 600115, Tamilnadu, INDIA

[www.helini.in](http://www.helini.in)

[info@helini.in](mailto:info@helini.in)

[helinibiomolecules@gmail.com](mailto:helinibiomolecules@gmail.com)

+91-44-24490433 / +91-9382810333